


1D NMR Acquisition Quick Guide



1 Getting Started:

- Logon using your Username and Password.
- Click on the VNMR Icon()
- Type *e < rtn >*. Place sample in spinner. Gauge properly. Place on top of magnet. Type *i < rtn >*.
- Click **Setup=>H1,CDCI3** (or desired *expt.*).
- Type *fixshims < rtn >*. Wait for beep.
- Click **acqi**. A new window will pop up.
- In new window, click **LOCK**.
- Check spinning. It should be 20.


4 Acquiring Your Spectrum:

- Type *nt=32* (or desired scans in multiples of 4) *< rtn >*. For ^{13}C , use *nt=1e6* *bs=8*.
- Type *ga < rtn >*. For ^{13}C , after a few data blocks are complete (message: BS # completed), type *wft* to process. When sufficient S/N is obtained, stop with *sa('bs')*.
- When complete, type *f full aph vsadj < rtn >*.

Manual Phasing (Optional):

- If autophasing did not work, type *lp=0 rp=0 < rtn >*.
- Click on **Phasing**. Using the *left* mouse button, click and hold on the **Rightmost** peak. Drag the mouse up or down to phase that peak.
- Using the *right* mouse button, click and hold on the **Leftmost** peak. Drag the mouse up or down to phase.
- Click **box**. If phasing is still poor, repeat above process.

2 Establishing Lock and Shimming:

- Click **LOCK off**.
- Move **lockpower** and **gain** slider to maximum.
- Move **Z0** slider slowly until one 'beat' is visible. 
- Click **LOCK on**.
- Reduce **lockpower** in units of 4 until correct value: 5-10 for D_2O , acetone- d_6 ; 10-15 for C_6D_6 ; 20-32 for CDCl_3 .
- Adjust **lockphase** to maximize lock level.
- Click **SHIM**.
- Adjust **Z1C** -1+ to maximize lock level. Repeat on **Z2C** -1+. Continue to alternate until no improvement. (Skip if no **Z1C**, **Z2C**).

3 Shimming (continued):

- Adjust **Z1** -64+ until maximum. Repeat with **Z2**.
- Adjust **Z1** -16+ until maximum. Repeat with **Z2**.
- Adjust **Z1** -4+ until maximum. Repeat with **Z2**.
- Click **CLOSE** to exit the lock/shimming window.
- Type *nt=1 ga < rtn >*. (Use with ^1H spectra only.)
- When complete, type *f full aph < rtn >*, expand around solvent peak or suitable well-resolved singlet. Type *vsadj < rtn >*.
- Is this peak well shimmed (i.e. is it narrow, symmetric, etc.)? If yes, proceed to acquisition. If not, click **acqi** then **SHIM**.
- Adjust appropriate shim (e.g. **Z1** for symmetric broadening or **Z2** for asymmetric peak shape).
- Repeat single scan acquisition. Reshim, if necessary.

5 Referencing Your Spectrum:

- Type *dscale < rtn >* and find your solvent peak (use the reference chart below or one near the instrument).
- With *left* mouse button click to the *left* of solvent peak.
- With *right* mouse button click to the *right* of solvent peak.
- Click **Expand**.
- Place cursor on top of solvent peak.
- Type *nl rl(< your solvent ppm > p) < rtn >*. For example, for CDCl_3 you would type *nl rl(7.24p) < rtn >*.

Common Deuterated Solvents:

CDCl_3	7.24p(^1H) 77p(^{13}C)	D_2O	4.63p(^1H)
Acetone- d_6	2.04p(^1H) 29.8p(^{13}C)	CD_3OD	3.30p(^1H) 49.0p(^{13}C)
Benzene- d_6	7.15p(^1H) 128.0p(^{13}C)	CD_2Cl_2	5.32p(^1H) 53.8p(^{13}C)
DMSO- d_6	2.49p(^1H) 39.5p(^{13}C)	DMF- d_7	2.91p(^1H) 35.2p(^{13}C)
CD_3CN	1.93p(^1H) 1.3p(^{13}C)		2.74p(^1H) 30.1p(^{13}C)

1D NMR Acquisition Quick Guide



6 Integrating Your Spectrum (Not for ¹³C NMR):

- Type *f full* <rt>. Click **Part Integral** or **Part Int**.
- Type *cdc dc cz* <rt>.
- Expand around first desired integral region.
- Click **Resets**.
- Use a **left** mouse click for each integral reset point. If you make a mistake, use the **right** mouse button to undo last reset point. To restart, type *cz* <rt>.
- Click **Full**, expand the next integral region, click **Resets**, left click your reset points, repeat for every region.
- When complete, click **Full**.

9 Plotting Your Spectrum:

- Typical example, *pl pscale ppf pir pltext page* <rt>.
- Type *ds* <rt>, expand desired plot regions, and repeat plot command.

Common Plotting Commands:

<i>pl</i>	plot spectrum
<i>pscale</i>	plot scale
<i>pir</i>	plot integral regions
<i>ppf</i>	plot peak frequencies
<i>pll</i>	plot line list with freqs in Hertz
<i>pltext</i>	plot text
<i>pltext(150,150)</i>	plot text in top right (<i>use with pll</i>)
<i>pap</i>	plot all parameters
<i>page</i>	send plot to printer

7 Referencing Your Integrals:

- Expand around Integral to be referenced.
- Place cursor on an integral region. The red vertical cursor must be on an integral trace.
- Click **Set Int** and enter an integral value (e.g. 2.0 <rt>).
- Type *vp=15* <rt>.
- Type *f full dpir* <rt> to display your integrals.
- Type *ds* <rt> to get an interactive spectrum.

Saving Your Data:

- Type *svf('your filename')* <rt>.

10 Logging Off of a Session:

- Type *e* to eject sample.
- Place standard in spinner. Gauge properly. Place on top of magnet.
- Type *i* <rt>.
- Type *exit* <rt>.
- Click the **EXIT** button on screen bottom, click **OK**.

8 Peak Picking:

- Type *f full ds* <rt>. Click **Th** and place yellow threshold line below top of smallest desired peak.
- Type *dpf* <rt>. If too many peaks, type *ds* <rt>, click **Th** and move threshold up. Type *dpf* <rt> to recheck.

Add Text to Your Spectrum:

- Type *text('your text here\\more text on new line\\even more text on a third line')* <rt>.

11 Analyze Your Spectrum:

